

Spe-ed™ Amberlite® XAD Resins



These analytical grade resins are prepared through a multi-solvent precleaning and rehydration. Each cartridge is slurry packed and the resins (cartridges or bulk) are shipped wet and ready to use.

Standard

Sorbent	100mg 1mL	200mg 3mL	500mg 3mL	500mg 6mL	1000mg 6mL	1000mg 12 mL
Amberlite XAD-4	2091	2092	2093	2096	2097	2099
Amberlite XAD-7	2191	2192	2193	2196	2197	2199
Amberlite XAD-16	2291	2292	2293	2296	2297	2299
Amberlite XAD-2000	3264	3265	3266	3267	3268	3269
Amberlite XAD-1180	3246	3247	3248	3249	3250	3251
Amberlite XAD-2010	3282	3283	3284	3285	3286	3287
Cartridges per box	100	50	50	30	30	20

Spe-ed™ LRC

(15mL reservoir cartridge)

	100mg	300mg	500mg
	3441	3443	3445
	3471	3473	3475
	3481	3483	3485
	3278	3279	3280
	3260	3261	3262
	3296	3297	3298
	30	30	30

SuperSpe-ed

Sorbent	2g 6mL	2g 12mL	2g 35mL	3g 12mL	3g 35mL	5g 35mL	10g 60mL	20g 60mL
Amberlite XAD-4	2098	5091	5092	5093	5094	5095	5097	5098
Amberlite XAD-7	2198	5191	5192	5193	5194	5195	5197	5198
Amberlite XAD-16	2298	5291	5292	5293	5294	5295	5297	5298
Amberlite XAD-2000	3270	3271	3272	3273	3274	3275	3276	3277
Amberlite XAD-1180	3252	3253	3254	3255	3256	3257	3258	3259
Amberlite XAD-2010	3288	3289	3290	3291	3292	3293	3294	3295
Cartridges per box	30	20	20	20	20	20	12	12

Bulk Material

Sorbent	100g
Amberlite XAD-4	3492
Amberlite XAD-7	3493
Amberlite XAD-16	3494

Sorbent	100g
Amberlite XAD-1180	3263
Amberlite XAD-2000	3281
Amberlite XAD-2010	3299

Amberlite® is a registered trademark of Rohm & Haas Co.



www.appliedseparations.com

930 Hamilton St. • Allentown, PA 18101
610.770.0900 • 610.740.5520 (fax)

spe14 12/03